

LUNA BELT



Luna Belt is a portable pressure and volume ventilator, that combines the latest technological innovations with a easy-to-use functionality.

Luna is designed for the hospital and homecare.

The ventilator offers a wide range of pressure support and volumemodes. It can be used with all conventional interfaces (mask, mouthpiece) and with the PBAir® interface, through the IAPV mode.

The IAPV mode (Intermittent abdominal pressure ventilation) allows an effective non invasive daily ventilation, avoiding tracheotomy and common issues (facial decubitus) of the mask ventilation.

APPLICATIONS

- **IAPV** (intermittent abdominal pressure ventilation)
- **IV** (invasive ventilation)
- **NIV** (non invasive ventilation)
- **MPV** (mouth piece ventilation)

INTERNAL LI-PO BATTERY, up to ten hours of therapy

FIELD OF APPLICATION

- Hospitals
- Home Care
- Respiratory Therapy Clinics

SUITABLE FOR

- Ideal for patients with acute/chronic respiratory failure
- Not designed to be used as an emergency ventilator
- Suitable for a wide range of patients' needs
- Provides continuous or intermittent ventilatory support for the care of subjects requiring mechanical ventilation
- Adult and pediatric up to 5kg

OUTSTANDING FEATURES

- Wide therapeutic versatility: Invasive and non-invasive ventilation, with active valve circuits and leakage circuits
- Accurate volume and trigger sensitivity
- Wide range of settings and alarms for each specific requirement
- "Easy to use" interface
- Small sizes, light weight and 10 hours of battery autonomy
- PBAir® interface, through IAPV mode